

Introduction

- "β-eye", is a dedicated planar coincidence camera suitable for *in vivo* molecular imaging of biomolecules and nanoparticles.
- β-eye is a unique benchtop system for whole-body mouse imaging. Its $5 \times 10 \text{ cm}^2$ field-of-view (FOV) allows static and fast dynamic studies.
- β-eye is the only truly portable planar coincidence system, offered in a safe suitcase with all components and ready for immediate use.
- β-eye fulfil the gap between *ex vivo* biodistributions and advanced multimodal imaging systems.
- Have a full biodistribution dataset, for all time points post injection, non-invasively, using only one animal.
- Obtain whole body images easily, right from the first second post injection.
- Speed up *ex-vivo* biodistributions by imaging all organs in a single view.
- Provide semi-quantitative planar images of significantly higher resolution, with no penetration depth limit at a much lower cost compared to optical imaging.
- Optimize imaging protocols.




System's Characteristics

Hardware

- "β-eye" camera
 - Dimensions: $40(L) \times 35(W) \times 30(H) \text{ cm}^3$
- Scintillation detection: pixellated BGO
 - FOV: $5 \times 10 \text{ cm}^2$
 - 22×44 elements
 - $2 \times 2 \times 5 \text{ mm}^3$ pixel size and septa 0.25 mm
- Photodetection: H12700A PSPMTs (Hamamatsu, Japan)
 - 1×2 arrangement - $5 \times 10 \text{ cm}^2$ FOV
- Number of detectors: 2
 - Separation distance: 6 cm
- FPGA electronics for data acquisition and processing
- Standard laptop

Software

- Fully comprehensive, user-friendly software
 - Real-time imaging
 - Database archive
 - Post-processing analysis
 - Reporting tool
 - DICOM export format



Packaging The system is delivered in a portable suitcase where all components are stored (mouse beds, phantoms, cables, laptop, power supply). The suitcase is safe for transportation by all means (airplane, bus, train) considered as standard luggage.





System's Performance

Table 1: Performance parameters


PARAMETER	350-700keV energy window
Spatial resolution	1.1 mm (CFOV)
Energy resolution	17 % @511 keV
System sensitivity	14 kcps/MBq
Timing resolution	2.2 nsec

Capillary (60 mm long, 0.6 mm diameter) filled with FDG 2 capillaries filled with FDG; Separation distance: 3 mm




Animal imaging

Mouse injected with 30 μCi FDG; 10 min post injection image



Mouse injected with 30 μCi FDG; 1 hour post injection image



Why "β-eye" is the right choice

Technology

- Low-cost benchtop system
- Easy versatile transportation
- Robust technology
- Semi-quantitative information
- Long-term operational system
- No special room requirements
- No need for technical staff
- User friendly software

Applications

- Whole-body dynamic studies
- Fast screening of promising biomolecules before detailed studies
- Dynamic studies for determining best biodistribution time-points
- Quality control imaging prior to *ex vivo* biodistributions
- Quality control pre-screening before multimodal imaging

Acknowledgement

This research has been co-financed by the Greek national funds through the National Strategic Reference Framework (NSRF) - 2014-2020, code: T1EΔK-01159 with acronym NAVIGATE.