



BIOEMTECH

Embracing scientists translate ideas into outcomes



γ -eye

A highly sensitive, benchtop, *in vivo* imaging system for all SPECT isotopes

General

Our vision is to accelerate preclinical research, towards clinical translation for promising drugs, through our high-quality services and products.

We are a strong research partner who guides, consults, and supports all preclinical research studies of our collaborators.

- As a manufacturing company of novel breakthrough imaging systems (*eyes*), we offer simplicity, speed, and efficiency on a daily workflow, during the first steps of testing novel compounds.
- As a preclinical CRO, we offer a one-stop-shop at our state-of-the-art Laboratories that covers a full chain of preclinical studies, following a Good Laboratory Practice approach in the daily routine.

BIOEMTECH's ultimate goal is to aid scientists innovate, at every step of their research.




BIOEMTECH

L. Mesogeion 387, Athens, Greece


info@bioemtech.com • www.bioemtech.com

Introduction

BIOEMTECH's γ -eye™ is a novel imaging concept, developed specifically for accelerating the early stages of preclinical research. Suitable for a wide range of SPECT isotopes, γ -eye™ enables real time *in-vivo* imaging of mice from time zero post-injection while also it can provide accurately all temporal and spatial physiological distribution characteristics of the studied compounds. Designed based on the end-user's needs and integrated in an easy-to-use and intuitive environment, γ -eye™ is a unique imaging tool that combines high flexibility, efficiency, and accuracy for an extensive range of applications.

With a footprint of just 44 cm × 46 cm x 40 cm and a weight lower than 40 kg, γ -eye™ is a truly desktop device that can turn any space into an imaging lab. γ -eye™ comes with a laptop PC, which serves for data acquisition and processing. Standard licenses of the complete software suite Visual | eyes, are included.

Technology – Specifications



A. General Information

γ -eye™ is a novel screening tool that offers the prosecution of fast, efficient and with high precision imaging studies. A large variety of SPECT isotopes as well as other imaginable signals resulting from alpha-emitting radionuclides can be studied in γ -eye™ within seconds, without the cost of time-consuming post-processing routines. γ -eye™ is a complete imaging solution, designed specifically to accelerate all stages of preclinical research.

| | |
|---------------------------|------------------------|
| Modality | SPECT |
| Anatomical mapping | Artificial X-ray |
| Active FOV | 50 mm × 100 mm |
| Photodetectors | Photomultipliers Tubes |
| Scintillators | CsI:Na |


B. Performance

γ -eye™, among other characteristics, infers good spatial resolution of 1.9 mm, energy resolution below 19%, while its high sensitivity enables real-time dynamic imaging with timeframes of 10 sec or lower (depending on the injection activity). Characterized by its high flexibility and based on the end user's needs, γ -eye™ can be specifically optimized for certain applications and imaging studies.

| | |
|---------------------------|---------------------|
| Time frames | Down to 1 sec* |
| Sensitivity | 341 cps/MBq |
| Spatial resolution | Up to 1.9 mm @ 0 mm |
| Energy resolution | Below 19% |
| Dynamic range | 30 keV – 500 keV |

C. Exchangeable collimators

γ -eye™, is designed to provide maximum flexibility to the user. Based on the type of the study, users can easily exchange collimators, choosing between a general purpose, a high-resolution, a high-sensitivity and a high-energy design. Specific designs can be provided upon request.




D. Anatomical Mapping

γ -eye™ integrates an advance Artificial Intelligence Algorithm designed to synthesize morphological X-ray images by translating standard photographic images of mice. Artificially produced X-ray mouse images can be superimposed with functional radioisotope 2D images to enhance overall anatomical information.



Optical Photo



Real X-ray




Artificial X-ray

E. Animal Handling

To preserve animal's welfare and health, γ -eyeTM employs standard inputs for gas anaesthesia, fully compatible with third party systems. Anaesthesia then is provided into the mask of the animals, all throughout the imaging study. In addition, the system infers a heated imaging stage, thus maintaining the temperature of the animal at the desired level. Upon request, vital signs of the animal can be monitored including heart rate, respiratory rate, body temperature and oxygen saturation.



| | |
|-------------------------------|--|
| Anaesthesia | Standard inputs for gas anaesthesia; compatible with third party systems |
| Heating | Heated stage for optimum body temperature |
| Vital signs monitoring | Upon request |



F. Footprint and connectivity

γ -eyeTM's footprint and standard digital interface connectivity can turn any space into an imaging lab. In addition, γ -eyeTM is characterized by simple power requirements and anaesthesia connections, allowing real time imaging and quantification inside a clean room, overcoming limitations, and facilitating scientists get great results under challenging conditions.

| | |
|-------------------------|-----------------------------------|
| Outer dimensions | 44 cm (L) × 46 cm (W) × 40 cm (H) |
| Weight | 40 kg |
| AC input range | 100-240 VAC |
| PC Connectivity | USB 2.0 Type A and GB Ethernet |
| Outer shielding | Sheet metal and acrylic |




G. Visual Eyes Software


VISUAL | eyes is a complete software suite that serves for image acquisition, system control and analytical post-processing routines. Among other possibilities, users through VISUAL | eyes software, can generate imaging studies using custom and/or pre-defined protocols, obtain quantitative information in user's defined Region of Interests and export images to DICOM format.

| | |
|---------------------------|---|
| Fast acquisition | Simplified procedure in a robust environment- |
| | Real time image visualization during the scan |
| Database | Raw data, DICOM storage, Compatibility with third party software |
| Anatomical mapping | Fusion with X-ray images artificially generated based on the mouse structural characteristics |
| Imaging protocols | Pre-defined and user's defined imaging protocols |
| Post processing | Integrated ROI manager for detailed post processing image analysis |
| License | Standard license for Mac and Windows |

Live imaging console



Post Processing suite



Indicative Studies

- **Tumor targeting for radiotherapy based on Pb212 radiopharmaceutical**




Figure 1. 14 uCi of Pb212 – 20min imaging time

- **Tumor targeting using I-125 radiopharmaceutical**




Figure 2. 40 uCi of I-125 – 1h dynamic scan

- **Inflammation imaging using Tc-99m compound**




Figure 3. 20 uCi of Tc-99m paw injection – 20min imaging time




BIOEMTECH

L. Mesogeion 387, Athens, Greece

info@bioemtech.com • www.bioemtech.com